





Lingine On Actation

Lab Simulation and Correlation to Engine Testing

Ricardo Hein
Conexo Inc.
Atlanta, GA
+1 678-792-4692
rhein@oilpas.com

Aeration properties of lubricating oils


Terms: Air entrainment, Microfoaming, dispersed gas

- DEFINITION: Aeration is a characteristic of an oil for its tendency to retain free air above the absorption limit at a given pressure and temperature.
- Oil aeration is important for good mechanical operation.
 Free air impairs mechanical performance because of:
 - 1. Cavitation in components by absorption and release
 - 2. Wear by irregular (potholed) oil film
 - 3. Oil compressibility ("sponginess") affects response in hydraulically driven component operation
 - 4. Oxidation reactions in a larger area of contact with oxygen and pressure cycles affects drain cycle and deposit/varnish formation

Aeration, Air Release, Foaming

- Previous test methods for aeration testing:
 - ASTM D-8047 Engine oil aeration resistance in CAT C13 heavy-duty diesel engine. COAT replaces D6894 as part of the PC-11
 - Air release test ASTM D-3427 and foaming test ASTM D-892 are insufficient to characterize air/oil interaction in performance drivetrain components
 - ISO 12152 Flender foaming test, stops test to measure increased volume in a ruler

Diagram of Testing Device




New test method to measure Oil Aeration

- Measures dynamic gas holdup in lubricating oils
- Principle of operation/test description
 - Aeration by simulation of a air vortex into the oil
- Aeration measured in steady agitation
 - 700 ml test fluid; Sample heating (25°C 120°C temperature control), time
 30 minutes
 - Turbine mixer develops vortex action to entrain air, Peristaltic pump moves a slipstream to test cell
 - OILPAS Optical imaging device measures bubbles in the test cell
 - The OILPAS algorithm evaluates online the dispersed gas ratio % in oil
 - Deaeration measured in time to release and speed of release
- Repeatability demonstrated

Set-up – Bench Test 2

Calculation the aeration % based on bubble measurement




⊙IL?∧**S**

Pump/Test Cell/ Microscope/Ima ging

Aerator, conditions control Speed, Temperature, Gas can be simulated in a dynamic cycle

Fluid Aeration Live Testing




Aerated Sample

Samples can be easily changed, allowing for scanning a large number of samples and formulations


Speed and Temperature can be adjusted to simulate a dynamic driving cycle

Adjustable frame for different sample sizes up to 5 liter


Software – Evaluation



Software – Data Diagrams




Aeration Study: Air Entrainment and Deaeration



Aeration Comparison 10 Engine Oils

All Samples are engine oil API SN grade 5W-30



Correlation to Engine Testing

Oil Sample A: 0W-20	Engine Test WOT	Aeration Bench	Comments
Mean Aeration During Test (%)	10%	9.1%	~ 10% correlation margin
Time for Deaeration (minutes)	Not reported	3:05	Time to release 99.8% of air

Oil Sample B: 5W-40	Engine Test WOT	Aeration Bench	Comments
Mean Aeration During Test (%)	19%	20%	~ 5% correlation margin !!!
Time for Deaeration (minutes)	Not reported	2:45	Time to release 99.8% of air

A direct correlation could be found between the bench test rig and the engine test at WOT regime

Conclusions

- Aeration screening in the lab can quickly benchmark the natural characteristic of an oil to retain free air
- Method was verified at testing a large number of different lubricating oils, additives and other fluids
- Measures both dynamic gas entrained and air release time,
- It shows correlation with ASTM D3427 test method
- Provides a measure of dynamic gas holdup and bubble size distribution (not provided by Coriolis sensor), thus adding further insight into the aeration characteristics of an oil.
- Provides OEMs/Equipment Builders and End Users an alternative method to test candidate fluids for aeration properties
- Test Method Guides OEM to design components considering the aeration characteristic of for the oil that is recommended for the specific application



Thank you for your attention. Any Questions?



2320 Starr Lake Drive

Acworth, GA 30101

USA

678 792 4692